


STANDING SLING WITH LOOP STRAPS

Instructions For Use



Device Overview:

DEVICE NAME:	DESCRIPTION:	PART NUMBER:	MATERIAL:	SWL:	
Standing Sling	Standing Sling with loop Straps Attachment	52460	Medium	Cotton/Poly Blend	250 kg. / 550 lbs.
		52462	Large		
		52464	Extra Large		



Washing Instructions:



WARNING!

Always read the instructions and the instructions for the lift being used with the sling.

Definition in this IFU:

Warning: Safety warning. Failure to understand or obey this warning may result in injury to you or to others.

Note: Important information for the correct use of this system or equipment.

Safety Instructions:

- Always check the label is readable and that the sling is in good condition without signs of wear and tear.
- Always make sure not to exceed the Safe Working Load (SWL) defined for the product.
- Always read the instruction manual for the lifting equipment being used with the sling.

Safe Working Load(SWL):

250 kg. / 550 lbs.

Intended Use:

The Standing Sling is a standard raising sling with anti-slip back and fixation belt providing great comfort when raising the user.

The Standing Sling is designed to be used with stand aid lifts and will assist the caregiver for raising and transferring a patient into a upright and/or sitting position. Loop Straps version is used with lifts/hoists and sling bars with two point loop attachment hooks.

Expected service life & Disposal:

Expected Service Life varies between 1–5 years, depending on usage and number of washing cycles. When washed once a week the sling has an Expected Service Life of 2 years. Expected Service Life is based upon a correct usage and washings procedure in accordance with this IFU.

Dispose of the sling according to your local regulations.

Before and After Every Use:

1. Always perform a risk assessment before starting the specific lifting session in order to make sure the product model is correct and that the lifting procedure is fully understood.
2. Before using the sling, read this instruction manual and the instruction manual for the lifting equipment being used.
3. When soiled or stained, the sling must be washed before use.
4. **Warning:** Check the sling for any fraying, loose stitching, attachments or any deteriorations.
5. After use, check the sling according to steps 3 and 4 and store in a dry area.

Service & Maintenance:

The sling must be checked by authorized personnel every 6 months according to ISO10535

Storing:

The sling should be stored in a dry and clean environment. Store in temperatures between (+32°-104°F)

Product care & Maintenance:

- Always follow the instructions on the product label and check the sling for deterioration.

Attachments:

The Standing Loop Sling will fit most stand aid lifts with standard loop attachment hooks. Always check lift instruction manual for details.

Customer Support/Trouble

Shooting:

- Please refer to the manual for the lifting equipment being used.
- If you have questions please contact Wy' East Medical about the product. (503) 657-3101



52602-R00

Standing Sling with loop straps is a raising sling with anti-slip back and fixation belt. It is designed to be used with sit-to-stand lifts and will assist the caregiver for raising and transferring a patient into an upright and/or seated position.

Applying the Standing Sling

1. Lean the patient slightly forward and place the sling with the anti-slip back against the patient and the label at the top. The bottom of the sling should be positioned at the patient's tailbone.
2. Connect the safety belt and pull it taut around the patient's waist.
3. Ensure both sides of the sling are positioned evenly around the patient.
4. Place the loop straps by feeding them under the sit-to-stand handlebars and secure them on the hooks. This will ensure the sling straps do not ride up into the patient's armpits.
5. Carefully perform the lift, cueing the patient to lean back into the sling.
6. Transfer the patient to the desired location.
7. Cue the patient to lean back into the sling before lowering.
8. Lower the patient until the sling straps are loose and remove them from the lift.