


WALKING SLING WITH LOOP STRAPS

Instructions For Use



Device Overview:

DEVICE NAME:	DESCRIPTION:	PART NUMBER:	MATERIAL:	SWL:	
Walking Sling	Walking Sling with Loop Straps Attachment	52466	Medium	Cotton/Poly Blend	180 kg. / 396 lbs.
		52468	Large		
		52470	Extra Large		



WARNING!

Always read the instructions and the instructions for the lift being used with the sling.

Definition in this IFU:

Warning: Safety warning. Failure to understand or obey this warning may result in injury to you or to others.

Note: Important information for the correct use of this system or equipment.

Safety Instructions:

- Always check the label is readable and that the sling is in good condition without signs of wear and tear.
- Always make sure not to exceed the Safe Working Load (SWL) defined for the product.
- Always read the instruction manual for the lifting equipment being used with the sling.

Safe Working Load(SWL):

180 kg. / 396 lbs.

Intended Use:

The Walking sling is designed with maximum comfort and safety in mind for the end user with some body control during exercise and rehabilitation purposes and provides easy installation for the caregiver and supportive weight relief for the user.

The Walking sling can be used with most mobile lifts/hoists, gait training and weight relief systems with Sling bars and loop attachment hooks.

The Walking sling is designed for use with lifts/hoists and Sling bars with two or four point attachment hooks.

Note: The Walking Sling is designed for support during training and not for lifting.

Expected service life & Disposal:

Expected Service Life varies between 1–5 years, depending on usage and number of washing cycles. When washed once a week the sling has an Expected Service Life of 2 years. Expected Service Life is based upon a correct usage and washings procedure in accordance with this IFU.

Dispose of the sling according to your local regulations.

Before and After Every Use:

1. Always perform a risk assessment before starting the specific lifting session in order to make sure the product model is correct and that the lifting procedure is fully understood.
2. Before using the sling, read this instruction manual and the instruction manual for the lifting equipment being used.
3. When soiled or stained the sling must be washed before use.
4. **Warning:** Check the sling for any fraying, loose stitching, attachments or any deteriorations.
5. After use, check sling according to steps 3 and 4 above and store in a dry area.

Service & Maintenance:

The sling must be checked by authorized personnel every 6 months according to ISO10535

Storing:

The sling should be stored in a dry and clean environment. Store in temperatures between (+32°-104°F)

Product care & Maintenance:

- Always follow the instructions on the product label and check the sling for deterioration.

Attachments:

The Walking Sling will fit most lifts with standard Sling bars with loop attachment hooks. Always check lift instruction manual for details.

Customer Support/Trouble

Shooting:

- Please refer to the manual for the lifting equipment used.
- If you have questions please contact Wy' East Medical about the product. (503) 657-3101



52603-R00

The Walking Sling is designed for walking patients with a dependent floor-based or overhead lift with a loop-style sling bar.

Applying the Walking Sling

1. Place the sling behind the patient with the label on the outside and the sling's center line aligned with the patient's spine.
2. Connect the sling belt's twin buckles and pull taut.
3. Position the leg straps under the patient's thighs. Connect each thigh buckle and pull tight.
4. Attach the sling straps to the sling bar.
5. Raise the sling bar to desired position and weight relief level.
6. Slowly cue patient to walk, checking that the sling straps stay attached and patient is comfortable.
7. After walking exercise has ended, lower sling bar until all straps are loose.
8. Remove loop straps from the sling bar.