

## PATIENT SLING COMPATIBILITY STATEMENT

### Information for Using Wy'East Medical Slings with other Manufacturer's Lift Equipment


Wy'East Medical Corporation has designed and manufactured safe patient handling products since its inception in 1987. Our company is registered with the Food and Drug Administration and is subject to periodic audits by the FDA to verify our compliance with Good Manufacturing Practices. Our slings are sold through national distributors and manufacturers of patient lifting equipment and are routinely used on lift equipment that is not manufactured by Wy'East Medical.

### Lift-Sling Compatibility


Wy'East Medical lift slings are fully compatible with other manufacturers' floor and ceiling lift systems using hook style sling bars. Lifting slings may be used with 2 or 4 point hook style sling bars. Repositioning slings may be used with multiple point hook style sling bars suitable for the size of sling being used. Wy'East Medical floor lifts are fully compatible with other manufacturers' slings utilizing the sling loop attachment method of attaching slings to floor lifts and ceiling lifts. The care provider must follow the responsible vendor's instruction manuals and recommendations with respect to the use, care, cleaning and inspection of slings and patient lifts.



Loop style Attachment



Loop attachment on hook style-  
loop sling compatible sling bar.



Clip style attachment -  
NOT COMAPATIBLE

Wy'East Medical recommends adherence to the **AASPHM Healthcare Recipient Sling and Lift Hanger Bar Compatibility Guidelines** for proper sling and lift inspection and safe consumer use; available at <http://aasphm.org/wp-content/uploads/AASPHM-Sling-Hanger-Bar-Guidelines-2016.pdf>

### Sling Liability

Wy'East Medical stands behind the safe application of our loop slings and accepts responsibility if a product defect results in an error in the practical application of the Wy'East Medical sling.

Wy'East Medical maintains product liability insurance on all the products we manufacture. Wy'East Medical will not be responsible for an error resulting from the improper use or application of the sling and/or lift.

### Sling Warranty

Wy'East Medical Corporation warrants to the original user that the patient sling manufactured by the company is free from defects in material and workmanship under normal use and service. The term of warranty from the date of purchase/shipment is two (2) years. During this term Wy'East will provide to the user a replacement patient sling.

In order to initiate a warranty claim the original user must provide prompt notice to Wy'East Medical Corporation within the applicable warranty period specifying the nature of the claim. Wy'East will issue instructions as to the procedure to follow. Wy'East reserves the right to confirm the existence of any claimed defect and to determine whether any damage was caused by neglect, misuse, unauthorized alteration or failure to follow maintenance instructions. With respect to a claim covered by this warranty, the defective part will be replaced or repaired at no charge, provided the defective part is returned to Wy'East Medical. The user has the responsibility for shipping charges and installation. Assistance is available by calling (503) 657-3101 in Oregon. This is the exclusive remedy for a claim under this warranty and in no event shall Wy'East be liable for consequential damages. THIS WARRANTY IS IN LIEU OF ALL OTHER WARRANTIES EXPRESSED OR IMPLIED.