

# AERIAL600

Maximum Assist Mobile Lift

WY'EAST  
MEDICAL



The Aerial is a maximum assist mobile patient lift. This diverse lift can be used for many patient handling tasks in addition to transfers, including floor recovery and physical therapy. The Aerial is light-weight, has 2 sling bar sizes and a patient scale available. Its low base option is convenient for those facilities with power-driven stretchers and beds with low under bed clearances.



**Patient Transfer**



**Edge of Bed**



**Therapy**

## Specifications

- Maximum Load 600 lb
- Number of Lifts per Charge 35-40
- Turning Diameter 64.3"
- Lifting Speed Unloaded 1.4 in/s
- Leg Under bed/Stretcher Clearance 2.5" or 3.8 "
- Standard Door Width Clearance 31.6"
- Standard or Small Sling Bar Options 24" or 15"



1-18001-R02



Patient Scale\*



2 sling bar sizes\*



Interchangeable  
Battery



QR Coded User  
Manual

\*optional accessories

## Features

Inter-sling compatibility

Low base option

Padded sling bar

Lightweight, aluminum frame

Encapsulated casters

800.255.3126 | wyeastmedical.com