




TECHNICAL SPECIFICATIONS:

PARAMETER	SPEC
• SWL 	500 lbs
• Device Part Number	51750
• Unit of Measure	1 Each
• Battery	Li-Ion 25.9V 2.25Ah
• Number of Lifts per Charge	35 - 40
• Charge Time	24 hrs fully discharged
• Total Weight	124 lbs
• Weight Loose Parts (Footrest)	12 lbs
• Wheel Diameter, Front	3"
• Wheel Diameter, Rear	4"
• Turning Diameter	46.1"
• Leg Under Bed/Stretcher Clearance	3.4"
• Standard Door Width Clearance	23.7"
• Footplate Dimensions	14" x 17" x 14"
• Lifting Speed Unloaded	1 in/s
• Frame	Steel & Aluminum
• FDA Device Class	Class 1
• FDA Product Code	FSA
• FDA Regulatory Number	880.5510



ENVIRONMENTAL CONDITIONS

	FOR SHIPMENT AND STORAGE	FOR USE
Temperature	14°F to 104°F (-10°C to 40°C)	41°F to 104°F (5°C to 40°C) ambient
Relative Humidity	20% to 80% non-condensing	20% to 80% non-condensing
Atmospheric Pressure	70 kPa to 106 kPa	70 kPa to 106 kPa

Note: Li-Ion batteries must be stored without exposure to direct sunlight



52701 R08 / AN 210617