

MATTRESSES

TotaLift / TC-300 Treatment Chair

WY'EAST
MEDICAL

Mattress Covers



51519

- Anti-static
- Anti-microbial
- Latex-free
- Fluid resistant
- 2.5" x 23" x 70"
- 9.4 lbs

COMFORT MEDICAL FABRIC: This durable pink cover is made of a patented thermoplastic fabric. Provides comfort and support for patients. Can be cleaned with most hospital-grade disinfectants.




51524

- Anti-microbial
- Breathable
- Latex-free
- Fluid resistant
- 2.75" x 23" x 70"
- 9 lbs

PRESSURE MANAGEMENT FABRIC: This cast-coated polyurethane polyester fabric stretches to conform to the patient and lowers interface pressure. Provides pressure relief for patients at higher risk for fragile skin such as ICUs, Neurology, and Oncology. Can be cleaned with most hospital-grade disinfectants.

Mattress Foam



1-51559-R08



51519 comes standard on the TotalLift. (51524 is available as an upgrade on the TotalLift).



51524 comes standard on the TC-300.